



Elevate your Workflow

SRS



Intracranial, Head & Neck



Craniospinal



Breast



SBRT



Lung



Abdomen



Pelvis



positioning patients for life.®

www.Qfix.com



Alta™ Standard Package

RT-4800-100

Alta™ Silverman Adapter Plate

Alta™ ArmShuttle™ Elite

InfinityLock™ Hand Grips

MR-Safe Locating Bars

Alta™ Multipurpose Device

SofTouch™ Knee Wedge Small

Qfix® Wall Mount

SofTouch™ Knee Wedge Large

SofTouch™ Foot Block

Quickly adapt to various clinical setup needs for different treatment types

Alta™ SBRT Accessories

RT-4800-B1

InfinityLock™ SBRT Compression Bridge (small)

RT-4800-B2

InfinityLock™ SBRT Compression Bridge (large)

RT-4558KV4-01

Low-Profile Compression Paddle

RT-4558KV4

Compression Paddle



Alta™ SRS Accessories

RT-4600-03-MRI

MR SRS Adapter



E0 Headrest

E2 Headrest



E4 Headrest

(SRS Headrests included with SRS Adapters)

Simplicity and versatility working together

At just 5.4kg (12 lbs), the Alta™ Multipurpose device is lightweight and therapist-friendly. Using Alta™ with compatible accessories facilitates device consolidation for patient positioning during radiotherapy.

The Alta™ spine area is clear and free of any objects and features a homogeneous platform with 3.8 mm WET for precise dose delivery, including hypofractionated treatments.



InfinityEdge™ & InfinityLock™



Alta™ Multipurpose Device provides a versatile platform that features the InfinityEdge™ for quick attachment of InfinityLock™ positioning accessories to streamline radiotherapy MR and CT imaging and treatment workflows.

SRS Accessories




RT-4600-03
SRS Adapter




RT-4600-03-MRI
MR SRS Adapter



E0 Headrest E2 Headrest E4 Headrest
(SRS Headrests included with SRS Adapters)



THERMOPLASTIC – OPTION 1



RT-B889KYCF2
Encompass™ SRS Fibreplast® System
(5) Fibreplast® VariablePerf™ Open View Mask
(5) Posterior Thermoplastic
(5) Optional IntegraBite™

THERMOPLASTIC & CUSTOM CUSHION – OPTION 2



RT-B889KYCF2A
(5) Anterior Only Fibreplast® Mask
(5) MOLDCARE® U Head Cushion
(5) Optional IntegraBite™



RT-4600-20
SRS Cushion Support (for forming MOLDCARE® U Head Cushion)



SBRT Accessories



RT-4800-B1
InfinityLock™ SBRT
Compression Bridge (small)

RT-4800-B2
InfinityLock™ SBRT
Compression Bridge (large)

RT-4558KV4-01
Low-Profile Compression Paddle

RT-4558KV4
Compression Paddle



RT-4558CB02
ZiFix™: Abdominal Thoracic
Motion Control System



Qfix
Avondale, PA
+1 484.720.6053
sales@Qfix.com
www.Qfix.com

2008576_EN

This brochure contains information about products which may or may not be available in particular countries. Each country has specific laws and regulations governing the commercialization of medical devices and the communication of information regarding medical devices in printed and digital media. Products appearing in this brochure or on our website may or may not have received approval, clearance, or marketing authorization by a governmental regulatory body in any particular country, or may have received approval, clearance, or marketing authorization for different indications and restrictions in different countries, or may be for investigational use only. A product's appearance in this brochure should not be construed as a solicitation or promotion for said product, nor any indication which is not authorized by the laws and regulations of the country where the reader resides. © 2020 All rights reserved. Qfix and Fibreplast are registered trademarks of Qfix Systems, LLC. Alta, ArmShuttle, Encompass, InfinityEdge, InfinityLock, IntegraBite, SofTouch, VariablePerf, and ZiFix are trademarks of Qfix. MOLDCARE is a trademark of Alcare, Co, Ltd.